

## Chromatin Beads

Thio-dCTP (100 mM), 1  $\mu$ mol (N-8002-1 from TriLink BioTechnologies)  
 Thio-dGTP (100 mM), 1  $\mu$ mol (N-8003-1 from TriLink BioTechnologies)  
 Biotin-14-dATP (0.5mM), 50 nmol (19524-016 from Invitrogen)  
 Biotin-21-dUTP, 0.5 mM (635701 from BD Biosciences)

DNA Polymerase I, Large (Klenow) Fragment, 5 Units/ $\mu$ l (M0210S from NEB)

Dynabeads M-280 Streptavidin 10mg/ml ( $6-7 \times 10^8$  beads/ml)<sup>1</sup>

Dynal® kilobaseBINDER™ Kit (SKU# 601-01, Invitrogen)

### Linearize plasmid

With BamHI and NotI, ideal ~ 4kb

### DNA biotinylation

#### 1. Fill-in reaction

DNA	30 $\mu$ g
Klenow buffer 10x	7 $\mu$ l
0.4 mM Biotin-14-dATP	8.7 $\mu$ l
0.5 mM Biotin-21-dUTP	7 $\mu$ l
5 mM Thio-dCTP	1 $\mu$ l
5 mM Thio-dGTP	1 $\mu$ l
Klenow (5 $\mu$ / $\mu$ l)	4 $\mu$ l
ddH <sub>2</sub> O	add 70 $\mu$ l

1h @ 37°C

#### 2. Separate unincorporated nucleotides

Precipitate DNA with EtOH

Check DNA concentration by measuring OD<sub>260</sub>

### Coupling DNA to beads

using the Dynabeads kilobaseBINDER™ Kit  
 (>182  $\mu$ g of a 4 kb DNA fragment/mg beads)

1. Resuspend the Dynabeads M-280 Streptavidin by shaking the vial to obtain a homogenous suspension.
2. Transfer 60  $\mu$ l (3  $\mu$ l beads/  $\mu$ g DNA) resuspended Dynabeads to an Eppi, place the tube in a Dynal MPC for 1 min until the Dynabeads have settled on the wall.

3. Remove supernatant. Avoid touching the Dynabeads pellet!
4. Add 240  $\mu\text{l}$  binding solution<sup>2</sup> and gently resuspend the beads.
5. Place the tube in the Dynal MPC and remove supernatant.
6. Resuspend the pellet in 134  $\mu\text{l}$  binding solution.
7. Add 134  $\mu\text{l}$  (20  $\mu\text{g}$  DNA) biotinylated DNA to the beads, mix carefully to avoid foaming!
8. Incubate at RT (18°C) o/n on a roller (about 1 turn / 10 s).
9. Place the tube in the Dynal MPC and remove supernatant (measure OD<sub>260</sub>).
10. Wash twice with binding solution.
11. Repeat a long (>4h) incubation at 4°C with 40% of the initial amount of DNA in “fresh” incubation mix. Measure OD<sub>260</sub> of the supernatant.
12. Wash the Dynabeads/DNA-complex twice in washing solution<sup>3</sup> and once in CSF-XB.
13. Calculate binding of DNA.
14. Resuspend the Dynabeads/DNA-complex in CSF-XB to get final conc. of 1 $\mu\text{g}$  DNA/ 3  $\mu\text{l}$  beads.

### **Making chromatin beads**

1. Take 230  $\mu\text{l}$  CSF-extract containing 1  $\mu\text{M}$  Nocodazole.
2. Resuspend 10  $\mu\text{l}$  of DNA-beads (3 $\mu\text{l}$  beads / $\mu\text{g}$  DNA) in 150  $\mu\text{l}$  of CSF-extract containing 1  $\mu\text{M}$  Nocodazole. Incubate for 5 – 10 min at 20°C (rather 16°C).
3. Add Ca<sup>2+</sup> to final conc. of 0.6 mM.
4. Incubate 90 min at 20°C, rotating.
5. Re-arrest in M-phase by adding 80  $\mu\text{l}$  CSF-extract.
6. Aliquot IMMEDIATELY and freeze in LN.

20  $\mu\text{l}$  of chromatin beads should contain enough beads for resuspension in 60  $\mu\text{l}$  new CSF!

---

1) in PBS pH 7.4, containing 0.1% BSA, 0.02 NaN<sub>3</sub>

2)

3) 10 mM Tris-Hcl (pH 7.5), 1 mM EDTA, 2.0 M NaCl