

Preparing CSF-extract from *Xenopus* Eggs

This protocol is essentially as described by Murray in „*Cell cycle extracts*“, Methods Cell Biol. 1991;36:581-605.

Frogs

HCG, SIGMA (in ddH₂O, filter sterilize, keep at 4°C)

PMSG, Intervet (optional)

Syringes

Needles Ø 0.4 mm (e.g. Microlance™ sterile needles, Becton Dickinson)

Containers

Room / Incubator at 16°C

- 1st injection: 50 U of human chorionic gonadotropin (HCG, SIGMA)
- alternatively: 50-100 U of pregnant mare serum gonadotropin (PMSG, Intervet)
- 2nd injection: 500 - 1000 U of HCG, 16 hours prior to when eggs are needed
- keep frogs separated in containers in 1 x MMR at 16°C

Solutions

Stock solutions

Cytochalasin B (10 mg/ml in DMSO; use 1:1000) SIGMA C2743

Cycloheximide (100 mg/ml in DMSO/EtOH; use 1:1000) Sigma C7698

Complete EDTA-free Protease Inhibitor (in ddH₂O; use 1:100) Roche
(keep at -20°C)

EDTA 0.5 M

EGTA 0.5 M

MgCl₂ 1M stock

CaCl₂ 1M stock

(filter sterilize, keep at RT)

10 x MMR

1 M NaCl

20 mM KCl

10 mM MgCl₂

20 mM CaCl₂

1 mM EDTA

50 mM HEPES pH 7.8

store at RT

5 l

292.2 g

7.46 g

10.17 g

11 g

10 ml (0.5 M stock)

59.58 g

20 x XB-salts

2 M KCl

20 mM MgCl₂

1 l

149.12 g

20 ml (1M stock)

2 mM CaCl₂ 2 ml (1M stock)
 filter sterilize, keep at 4°C

Prepare freshly

<u>CSF-XB</u>	<u>1 l</u>
10 mM HEPES	2.38 g
50 mM Sucrose	17.11 g
100 mM KCl, 1 mM MgCl ₂ , 0.1 mM CaCl ₂	50 ml (20 x XB-salts)
5 mM EGTA	10 ml (0.5 M stock)
adjust to pH 7.7 with KOH	

<u>Dejelly Buffer</u>	<u>1 l</u>
2% (w/v) L-Cysteine	20 g
0.25 x MMR	25 ml (10 x MMR)
adjust to pH 7.8 with 5N NaOH	

Extract preparation

Plastic Pasteur pipettes (e.g. Peske # 18-4321; 3 ml Pasteur-Pipetten graduiert)
 Clinical centrifuge
 Beckman Ultra-Clear centrifuge tubes (13 x 51 mm; Reorder No. 344057)
 Sorvall HB-4 or HB6 rotor
 Polyallomer tube (Sarstedt)
 Syringes
 Needles Ø 1.2 mm

- Collect eggs and keep in separate batches if possible
- Wash eggs in 1 X MMR
- Dejelly 2 x 2.5 min with Dejelly Buffer
- Wash 2 x with XB-CSF
- Wash 1 x with XB-CSF + protease inhibitors
- Preload Beckmann tube with 1 ml XB-CSF + 1 µl Cytochalasin B
- pack eggs (place Beckmann tube into polyallomer tube, spin 500 rpm for 30 s, then 2000 rpm for 1.5 min at 14°C)
- spin (rotor HB4 / HB6 RT!) in cold (4°C) centrifuge for 20 min at 10 700 rpm
- collect cytoplasmic fraction with syringe (Ø 1.2 mm)
- add cytochalasin B, final conc. 10 µg/ml
- add protease inhibitors
- (spin extract for add. 60 min at 50 000 rpm in TLS55, collect supernatant as high speed extract)

Testing the Extract

labeled tubulin
 any kind of nucleating source (sperm nuclei, centrosomes, RanGTP, TPX2,...)

fix
coverslides
coverslips

Fix

0.3 ml formaldehyde 37%
0.6 ml glycerol 80%
1 ml 1x MMR
1 μ l Hoechst No. 33342 (10 mg/ml)

100 μ l aliquots
keep at -20°C

test 20 μ l extract for CSF-arrest by adding

- labeled tubulin at 1/50 – 1/200 (depending on labeling ration)
- nucleating source
- incubate the extract for 20 min at 16-20°C
- fix 3-4 μ l reaction by squash fixing with 3-4 μ l fix solution under a coverslip
- check for mitotic structures

Interphase Extract

in case you need interphase extract, add

- 0.6 mM CaCl₂
- Cycloheximide, final conc. 100 μ g/ml
- incubate for 40 min at 16-20°C
- check for interphasic structures as described above